e * BMC Pre-Lab

BMC Pre-Lab

INCLUDE THIS SHEET AS THE FIRST PAGE OF YOUR REPORT.

Sissi Wong (Xigue Jung Lm (Don's) Lee

Before the 1st Day of Lab

Pre-lab Discussion Questions

first page of your report. Without it you will lose grade points. You should be prepared to discuss at least the following before you come to labs tis your responsibility to discuss this lab with an instructor before your first day of your scheduled lab period. This signed sheet must be included as the

1. What are the masses of the various nanoparticles you will be observing in the lab? How many atoms are in a single particle? What is the uncertainty in

plan to observe and simulate them in Matiab. (You should choose at least one particle 1 µm or larger and one smaller.) a use your simulated data to calculate the diffusion coefficient, D in each case. Explain how you arrived at your answer. 2. Using the microscope, you will observe a minimum of two different size particles in at least four different viscosity solvents. Choose the conditions you these numbers? Data sheets for the nanoparticles are available on the BMC Reprint List.

Keep these scripts. When applying your data you can create artificial data sets on which to test your analysis techniques ♣/What is the uncertainty of your estimate of D? How does it vary with the number of simulated data points? Explain your strategy for making observations in the lab. What additional sources of error (these are significant) will come in to play? How will you account for them? 2006

Staff Signature Afficia Completed on before the first day of lab? (circle) Yes / No

Mid-lab Questions

On day 3 of this lab, you should have completed the following: Show them to an instructor and ask for a signature

Mid-lab Questions Part I

- 2. Using a slide with a combination of 10µm and 0.44µm polystyrene spheres, show how to set up Köhler illumination.

 3. How many nanometers per pixel are captured at 10x, 20x, and 40x?

 4. Draw diagram of darkfield illumination. Explain how it is possible to see 40 nm objects with visible light (400-750 nm wavelengths).

 5. Set up dark-field illumination.

Completed by the second day of lab? (circle) Yes

Mid-lab Questions Part II

By day three of this lab, you should have collected some particle tracks and made several movies. Show one of the particle tracks to an instructor. What value of D did you calculate from the track? How close is this to the theoretical value? You can do this either with the BMC application or with the Matlab scripts. Show and explain your averaging and centroiding code. How do they work? see next signature

pleted by the third day of lab? (circle) Yes

Please also fill out the Student Evaluation of Experiment.