



QIE Pre-Lab

INCLUDE THIS SHEET AS THE FIRST PAGE OF YOUR REPORT.

Before the 1st Day of Lab [1]

Student's Name: Jung Lin (Doris) Lee

Partner's Name: Xoyue Wang

Before the Lab

Before using the apparatus in this experiment, you must complete training in the safe use of lasers detailed on the [Laser Safety Training](#) [2] page. This includes readings, watching a video, taking a Laser Quiz, and filling out a Laser Training Certificate form. You must turn in the completed forms into the 111 Lab Staff before you start the experiment.

Pre-lab Questions

It is your responsibility to discuss this lab with an instructor before your first day of your scheduled lab period. This signed sheet must be included as the first page of your report. Without it you will lose grade points. You should be prepared to discuss at least the following before you come to lab:

1. What is the significance of an observation of a violation of a Bell inequality? What does it mean for the validity of Quantum Mechanics as well as for the character of possible alternative theories?
2. A half-wave plate is a slice of birefringent crystal that allows incident light whose polarization is perpendicular to some preferred axis (the optical axis) to pass through unchanged but delays by 180 degrees, i.e., changes the sign of, light whose polarization is parallel to the optical axis. Given this, what is the output polarization of a half-wave plate whose optical axis forms an angle θ with the polarization axes of the incident light? If you wish to change horizontally polarized light into light that is in an equal superposition of horizontal and vertical polarizations—that is, polarized at 45 degrees relative to the horizontal—at what angle should you set your waveplate?
3. Say along each entangled photon path, we have the following optical elements, listed in the order the photon will encounter them: a half-wave plate, a horizontal polarizer, and a detector.
 - a. Four standard measurements we use to see what state our light is in are $N(0,0)$, $N(0,90)$, $N(45,45)$, and $N(90,90)$, where $N(\alpha,\beta)$ is the number of coincidences with the two detector wave plates set to rotate photon polarizations by α and β . **Note that these are *not* the actual waveplate settings.** (See previous question.) Assuming a perfect Bell state and zero phase difference between the vertical and horizontal downconverted photons, i.e. $|\psi\rangle = (1/\sqrt{2})(|hh\rangle + |vv\rangle)$, what values do we expect to see at each of these polarizations, if we know the system is generating 1000 downconverted photons per second? Show the instructor code to calculate the expectation values for arbitrary angles, it will be useful later to analyse the data on the fly.
 - b. Of course, the real use of these measurements is when our downconverted photons aren't in the ideal entangled case. Let's look at two important cases of this. Suppose that our first half waveplate is slightly off from its correct value, and as a result our down-converted pairs are more likely to be horizontally polarized than vertically polarized. For example, instead of being in the state $|\psi\rangle = (1/\sqrt{2})(|hh\rangle + |vv\rangle)$, they are now in the state $|\psi\rangle = (\sqrt{3}/2)|hh\rangle + (1/2)|vv\rangle$. How do we expect each of these measurements to change?

again, but they are no longer entangled—that is, it is emitting light that is a 50–50 mix of horizontally polarized photon pairs and vertically polarized photon pairs. What would you now expect from each of these measurements? This is what we could observe if our experiment followed a local hidden variable theory.

- d. Now suppose that our setup is producing a pure entangled (Bell) state again, but that there is a phase difference between the horizontal and vertical components of the state. Now what do we expect from each of these measurements? This is why we have an angled birefringent plate in the 405nm beam path, to compensate for the space between the two BBO crystals.

✓ 4. Our coincidence counting works in the following manner: each time a photon hits one detector, the FPGA waits 5ns afterwards to see if a photon hits any other detector. Of course, there are some background counts on each detector due to ambient light, and every now and then two of those counts will come close enough together to register a "false" coincidence. Determine the expected rate of these false coincidences in terms of the rate of counts on each detector and the coincidence time window. If each detector gets 10000 background counts per second, what rate of background coincidences do you expect for a 5-ns coincidence window? How does this scale if you double the amount of background light (which doubles the count rates on each detector)?

✓ 5. What is the statistical error on the number of counts on each detector? How about the statistical error on the number coincidences?

6. Suppose that as you start to perform your experiment, two of your friends come crashing into the room to see what you are doing, one through each wall (like the Kool-Aid Man) at a speed of .9999c (in opposite directions) when they enter. Unfazed, you explain the experiment and how it works, but once you start to describe how the wave function collapses when measured, one of your friends cuts you off.

"That's all fine," he says, "but there's just one thing: I saw you running those photon pairs as I arrived,"—he's quite a talented fellow—"and the photons going to detector A were always about 19 ns delayed from the ones going to detector B. So it was the measurement of those photons that collapsed the entangled state, right?"

Before you can answer, your other friend jumps in: "You mean the ones going to detector B were delayed, right?" she asks, puzzled. "That's what I saw."

✓ 7. How do you reconcile what your friends saw? What implications does this have for how entanglement works? Explain which practices you will follow while setting up and disassembling the optics as well as with the fibers (see Optics Tutorial [3]). In particular, show the instructor how you will plug a fiber into or unplug from a connector.

Staff Signature: _____ Date: 3-4-2016

Completed before the first day of lab? (circle) Yes / No

Mid-lab Questions and Sign Off Sheet

By day 3 or 4 of this lab, you should have successfully aligned the laser and produced a well tuned Bell state.

1. Record and report your half-wave plate calibrations.
2. Create a graph that shows the purity of your Bell state and show it to a GSI for a signature. Also, show your quantified bell state (coefficients C_1 , C_2 and phase ϕ).
3. Which pairs of half-wave plate angles will you use for your measurement of S ? Why are these the best values to use for your entangled state? Make sure you have math to back you up.
4. Please measure the output power of the laser diode after the optical isolator and show and tell the professor or GSI. Document this number.

Staff Signature: _____ Date: 3/10/16

Completed on the fourth day of lab? (circle) Yes / No

INCLUDE THIS SHEET AS THE FIRST PAGE OF YOUR REPORT.

Please also fill out the Student Evaluation of Experiment [4].

Powered by Drupal